



## 28623 DUFTDOS-NT35-R

### FUNCTIONAL DESCRIPTION

Peristaltic pump for dosing fragrances into the steam line of steam baths, with integrated intensity control. The dosing capacity is set via the „pause and dosing time“ on the five-line LC display using 3 input keys. Operation only with steam generator switched on. Dosing capacity maximum 3 l/h.

Optionally, the dosing pump can be operated with a separately available contact sensor. This measures the temperature at the steam line. As soon as a certain, adjustable temperature between 35 °C and 95 °C is reached, the steam pump is activated.

### SCOPE OF DELIVERY

- Control NT35 with LC display built into housing
- Peristaltic pump SR10 approx. 3/h
- Dosing hose 1,6 x 1,6 Viton
- Dosing cassette SR10 SR10
- Power cable with euro plug
- Suction hose PE 4 x 1 - 1,5 m



### TECHNICAL DATA

Electrical connection	100-230 V/AC 50-60 Hz
Connected load	10 W
Back pressure max.	ca. 0,5 bar
Size (W/H/D)	80 x 180 x 130 mm
Weight	approx. 1,4 kg

### ACCESSORIES / OPTIONS

- Dosing set consisting of: Dosing valve VA 3/8“ 4 x 1 and 3 m dosing hose PTFE 4 x 1 - Part No. 12596
- Console with wide-neck bottle 1,5 l - Part No. 13892
- Console for 10 l canister - Part. No. 14391
- Push button, Alu Dibond - Part. No. 19011
- 1 x Bonnet 1-fold - Part. No. 21030-1
- Steam bath push button - Part No. 23561
- Sealing flange - Part No. 22116-1
- Mounting sleeve 40/42mm - Part No. 20925
- Contact sensor PT 1000 - Part No. 25652